


Product Datasheet

ARPC1B antibody (orb18897)

Description	Goat polyclonal antibody to ARPC1B
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ARP2/3 subunit 1B
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. WB: Approx 37kDa band observed in Human Liver lysates (calculated MW of 40.9kDa according to NP_005711). Recommended concentration: 2-4 µg/ml. A minor band of unknown identity was also consistently observed at 28kDa.
Protein Sequence	AYHSFLVEPISCH
Clonality	Polyclonal
MW	40.9
Entrez	10095
Dilution Range	ELISA: 1:32000, WB: 2-4 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...



Immunohistochemical staining of Human Th...