



Product Datasheet

XAB2 antibody (orb18898)

Description	Goat polyclonal antibody to XAB2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	XAB2 / HCNP
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: In paraffin embedded Human Brain (Cerebral Cortex) shows nuclear staining. Recommended concentration, 3-4µg/ml.WB: Approx 90kDa band observed in human fibroblast A431 lysates (calculated MW of 100kDa according to NP_064581.2). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	SVPAAVFGSLKED
Clonality	Polyclonal
MW	100
Entrez	56949
Dilution Range	ELISA: 1:8000, WB: 0.3-1 µg/ml, IHC-P: 3-4µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis
of A431 lysate (RI...)



Immunohistochemical
analysis of formalin...