


Product Datasheet

AKAP8 antibody (orb18902)

Description	Goat polyclonal antibody to AKAP8
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	AKAP8 / AKAP95
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Experiments gave a band at approx 100kDa in Human Liver lysate after 1 µg/ml antibody staining. This band corresponds to the 95kDa band reported in literature (PMID: 10601332). This protein has a calculated MW of 76.1kDa according to NP_005849. Recommended concentration: 1-3 µg/ml.
Protein Sequence	DQGYGGYGAWWSAG
Clonality	Polyclonal
MW	76.1
Entrez	10270
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb18902
staining
(1ug/ml) of
Human Live...