

Product Datasheet

Swiprosin 1 antibody (orb18908)

Description

Goat polyclonal to Swiprosin 1

Species/Host

Goat

Reactivity

Bovine, Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Target

Swiprosin 1 / EFHD2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

This antibody does not cross-react with Mouse Swiprosin 2. Data kindly provided by Sebastian Duetting and Dirk Mielenz. ELISA: Peptide ELISA: antibody detection limit dilution 1:16000. IHC: In paraffin embedded human tonsil stains cells in the proximity of the epithelial cryptae, granulocyte neutrophils, macrophages within the germinal centre and in scattered cells (particularly in the interfollicular area and associated with epithelium) thaWB: Approx 28kDa band observed in Human Lymph Node and Tonsil lysates (calculated MW of 26.7kDa according to NP_077305.2). Recommended concentration: 0.1-0.3 µg/ml.

Protein Sequence

CAAFKELQSTFK

Clonality

Polyclonal

MW

26.7

Entrez

79180

Dilution Range

ELISA: 1:128000, WB: 0.5-1.5 µg/ml

Expiration Date

12 months from date of receipt.

250kDa
150kDa
100kDa

50kDa
37kDa
25kDa
20kDa
15kDa

Western blot analysis of human lymph nod...

250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
20kDa
15kDa

Western blot analysis of Human Thymus an...

250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

Western blot analysis of Human Thymus (L...