



Product Datasheet

DOK5 antibody (orb18915)

Description	Goat polyclonal to DOK5
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	DOK5
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.IHC: In paraffin embedded tonsil, stains scattered cells in the interfollicular area compatible with mast cells. Recommended concentration, 1-2ug/ml.WB: In transfected HEK293 transiently expressing DOK5 a band of approx. 38kDa is observed. This band is not observed in the non-transfected HEK293. The calculated molecular size is 35.5kDa according to NP_060901.2. Recommended concentration, 1-3 µg/ml.
Protein Sequence	CLHRTETFPAYRSEH
Clonality	Polyclonal
MW	35.5
Entrez	55816
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml, IHC-P: 1-2µg/ml
Expiration Date	12 months from date of receipt.



HEK293
overexpressing
DOK5
(RC211718)
an...



orb18915
(1ug/ml)
staining of
paraffin e...