


Product Datasheet

FOXN1 antibody (orb18923)

Description	Goat polyclonal antibody to FOXN1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	FOXN1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Peptide ELISA: antibody detection limit dilution 1:16000. IHC: In paraffin embedded Human Thymus shows nuclear staining of thymocytes in the densely populated cortex. Recommended concentration: 4-6µg/ml.
Protein Sequence	SVYLSPPSSKPVALA
Clonality	Polyclonal
MW	68.9
Entrez	8456
Dilution Range	ELISA: 1:16000, IHC-P: 4-6 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Th...