

Product Datasheet

ADH1A, B, C antibody (orb18937)

Description	Goat polyclonal antibody to ADH1A, B, C
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	pan ADH
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.IHC: Customer found this product to work in IHC on Human Liver.WB: Approx 38kDa band observed in Human Liver lysates (calculated MW of 39.9kDa according to NP_000658.1, NP_000659.2 and 000660.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	STAGKVMKCKA
Clonality	Polyclonal
MW	39.9
Entrez	124 , 125 , 126
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb18937
(1ug/ml)
staining of
Human
Live...