



Product Datasheet

AKR1C3 antibody (orb18943)

Description	Goat polyclonal antibody to AKR1C3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	AKR1C3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 35kDa band observed in human liver and human breast lysates (calculated MW of 36.9kDa according to NP_003730). Recommended concentration: 0.01-0.1 µg/ml.
Protein Sequence	CFASHPNYPYSDEY
Clonality	Polyclonal
MW	36.9
Entrez	8644
Dilution Range	ELISA: 1:32000, WB: 0.01-0.1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human breast ly...



Immunohistochemical staining of Human Sm...