

Product Datasheet

XAGE1 antibody (orb18946)

Description

Goat polyclonal to XAGE1

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications

Immunohistochemical
staining of Human
Te...

Target

XAGE1 (Isoform 1d) / GAGED2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

This antibody is expected to recognise only the short isoform represented by NP_597673.1.ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as small as the predicted size of approx. 8.0kDa according to NP_597673. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

CGFGFRRQGEDNT

Clonality

Polyclonal

MW

8

Entrez
9503
Dilution Range

ELISA: 1:8000, IHC-P: 3-6µg/ml

Expiration Date

12 months from date of receipt.