

Product Datasheet

FARP1 antibody (orb18947)

Description	Goat polyclonal antibody to FARP1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	CDEP / FARP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Preliminary experiments in Human Liver, Kidney and Lung lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	GEIEQRPTPGSRL
Clonality	Polyclonal
MW	119; 14.1; 122
Entrez	10160
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.