



Product Datasheet

LRRK2 antibody (orb18949)


Description	Goat polyclonal antibody to LRRK2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	LRRK2 / PARK8 (near C Terminus)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.IHC: In paraffin embedded Human Brain (hippocampus) shows staining of neuronal cytoplasm. Recommended concentration, 1-2ug/ml.WB: This product has been successfully used in both Western blot and immunohistochemistry on human brain: Alegre-Abarrategui J et al., Neuropathol Appl Neurobiol. 2007 Oct 26; PMID: 17971075. Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	CELAEKMRRTSV
Clonality	Polyclonal
MW	286
Entrez	120892
Dilution Range	ELISA: 1:128000, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Hi...



Immunohistochemical staining of Human Hi...



Immunohistochemical staining of Human Co...