

Product Datasheet

Centromere protein F antibody (orb18953)

Description

Goat polyclonal to Centromere protein F

Species/Host

Goat

Reactivity

Bovine, Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications
Target

Centromere protein F / mitotin

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the predicted size of approx. 358kDa according to NP_057427. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

GLESKGSENCVKVQ

Clonality

Polyclonal

MW

358

Entrez

1063

Dilution Range

ELISA: 1:64000, IHC-P: 3.75 µg/ml

Expiration Date

12 months from date of receipt.