

Product Datasheet

CENPE antibody (orb18954)

Description

Goat polyclonal antibody to CENPE

Species/Host

Goat

Reactivity

Bovine, Canine, Human, Rat

Conjugation

Unconjugated

Tested

ELISA

Applications
Target

Centromere protein E

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the predicted size of approx. 316kDa according to NP_001804. Therefore we cannot recommend an optimal concentration and the product is Investigative Grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

HASSGKDVPECKTQ

Clonality

Polyclonal

MW

316

Entrez

1062

Dilution Range

ELISA: 1:32000

Expiration Date

12 months from date of receipt.