


Product Datasheet

TRAK1 antibody (orb18962)

Description	Goat polyclonal antibody to TRAK1
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	OIP106 / TRAK1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 90kDa band observed in Human Heart and Human Lung lysates (calculated MW of 77.3kDa according to NP_055780). Recommended concentration: 0.3-1.5 µg/ml.
Protein Sequence	CGAKLSKQTSRLR
Clonality	Polyclonal
MW	77.3
Entrez	22906
Dilution Range	ELISA: 1:64000, WB: 0.3-1.5 µg/ml
Expiration Date	12 months from date of receipt.



orb18962
staining
(0.3ug/ml) of
Human He...