

Product Datasheet

BOLL antibody (orb18976)

Description

Goat polyclonal antibody to BOLL

Species/Host

Goat

Reactivity

Bovine, Canine, Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications

Immunohistochemical
staining of Human
Te...

Target

BOULE / BOLL

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. WB: Preliminary experiments gave an approx 22kDa band in Human Testis lysates after 1.5 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 32.6kDa according to NP_932074.1 and 31.3kDa according to NP_149019.1. The 22kDa band was successfully blocked by incubation with the immunizing peptide. However, literature reports a 36kDa band in Human and Mouse Testes (PMID: 11390979). Have any further splice variants/modified forms been reported

Protein Sequence

QPEPIKTVWSIHY

Clonality

Polyclonal

MW

32.6; 31.3

Entrez
66037
Dilution Range

ELISA: 1:32000, WB: 1.5 µg/ml, IHC-P: 4-6µg/ml

Expiration Date

12 months from date of receipt.

