


Product Datasheet

SRD5A1 antibody (orb18980)

Description	Goat polyclonal antibody to SRD5A1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	SRD5A1 / 5-alpha reductase 1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.
Protein Sequence	ATATGVAEERLLC
Clonality	Polyclonal
MW	29.5
Entrez	6715
Dilution Range	ELISA: 1:128000, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Prostate tissue showing SRD5A1 expression.