



Product Datasheet

NPM1 antibody (orb18984)


Description	Goat polyclonal antibody to NPM1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Target	Nucleophosmin / NPM1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 37kDa band observed in Human HeLa lysates (calculated MW of 32.6kDa according to NP_002511.1 and 29.5kDa according to NP_954654.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	QEAIQDLWQWRKSL
Clonality	Polyclonal
MW	32.6; 29.5; 25.3; 18.9
Entrez	4869
Dilution Range	ELISA: 1:8000, WB: 0.01-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Spleen ly...



Western blot analysis of NIH3T3 (lane 1)...



Immunofluorescence analysis of HeLa cell...