


Product Datasheet

RBP1 antibody (orb18987)


Description	Goat polyclonal antibody to RBP1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	RBP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.IHC: In paraffin embedded Human Pancreas shows vesicular staining in the acinar cytoplasm but oriented towards the interlobular space.. Recommended concentration, 5-10ug/ml.WB: Approx 16kDa band observed in Murine NIH 3T3 lysates (calculated MW of 15.9kDa according to NP_002890.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	VEGVVCKQVFKKVQ
Clonality	Polyclonal
MW	22.3
Entrez	5947
Dilution Range	ELISA: 1:4000, WB: 0.5-3 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH3T3 lysate u...



Immunohistochemical staining of Human Pa...



Immunohistochemical staining of Human Li...