



Product Datasheet

DPP10 antibody (orb18988)


Description	Goat polyclonal antibody to DPP10
Species/Host	Goat
Reactivity	Bovine, Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Dipeptidylpeptidase 10
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 90kDa band observed in human bone marrow lysates (calculated MW of 90.0kDa according to NP_065919.2 and of 90.1 according to NP_001004360.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	DEGHNVSEKSKYHL
Clonality	Polyclonal
MW	90.9; 90.1
Entrez	57628
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human bone marrow...



Immunohistochemical staining of Human Co...



Immunohistochemical staining of Human Ce...