


Product Datasheet

ARRB2 antibody (orb18990)


Description	Goat polyclonal antibody to ARRB2
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	Arrestin beta 2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 40kDa band observed in Human Brain lysates (calculated MW of 46.1kDa according to NP_004304.1 and 44.4 according to NP_945355.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	HDHIPLRPQS
Clonality	Polyclonal
MW	46.1; 44.4; 48.0; 47.3; 45.6; 24.3
Entrez	409
Dilution Range	ELISA: 1:64000, WB: 0.5-1µg/ml, IF/ICC: 10µg/ml, FACS: 10ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 (A), U251 ...



Immunofluorescence analysis of paraforma...



Flow cytometric analysis of paraformalde...