



Product Datasheet

BID antibody (orb18994)

Description	Goat polyclonal antibody to BID
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Target	BID
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 26+24kDa bands observed in Human fibroblast A431 lysates (calculated MW of 26.8kDa according to NP_932070.1 and 22.0kDa according to NP_001187.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	TYVRSLARNGMD
Clonality	Polyclonal
MW	26.8; 22.0; 11.3
Entrez	637
Dilution Range	ELISA: 1:8000, WB: 0.1-0.3 µg/ml, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



orb18994 (0.1ug/ml)
staining of A431
lys...



Immunofluorescence
analysis of MCF7
cell...