

Product Datasheet

RASA3 antibody (orb18999)

Description	Goat polyclonal antibody to RASA3
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	RASA3 / GAP1IP4BP
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Preliminary experiments in Human Brain lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). - have any results been reported with other antibodies/lysates
Protein Sequence	QYKRDKFKKTKYGSQ
Clonality	Polyclonal
MW	95.7
Entrez	22821
Dilution Range	ELISA: 1:32000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.