

Product Datasheet

SNAP25 antibody (orb19005)

Description

Goat polyclonal to SNAP25

Species/Host

Goat

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Target

SNAP25

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Brain (Cerebral Cortex) shows axon/dendrite staining leaving cell bodies relatively unstained. Recommended concentration, 2-4ug/ml.WB: Approx 22kDa band observed in Human Brain and Mouse Brain lysates (calculated MW of 23.3kDa according to NP_003072.2 and NP_570824.1). Recommended concentration: 0.01-0.03 µg/ml.

Protein Sequence

DEANQRATKMLGSG

Clonality

Polyclonal

MW

23.3

Entrez


6616

Dilution Range


ELISA: 1:128000, WB: 0.03-0.1 µg/ml, IHC-P: 1-2 µg/ml

Expiration Date


12 months from date of receipt.



Western blot analysis of Mouse Brain lys...



Immunohistochemical staining of Human Ce...



Western blot analysis of Human Brain (Ce...