


Product Datasheet

HTR2C antibody (orb19008)

Description	Goat polyclonal antibody to HTR2C
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Serotonin receptor 2C / HTR2C
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 50kDa band observed in lysates of Human EBV immortalised Lymphoblastoids, T cell MOLT4 and monocytic cell line U937 (calculated MW of 51.8kDa according to NP_000859.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	QVENLELPVN
Clonality	Polyclonal
MW	51.8
Entrez	3358
Dilution Range	ELISA: 1:64000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Lymphoblastoid

...