


Product Datasheet

CD2BP2 antibody (orb19020)

Description	Goat polyclonal antibody to CD2BP2.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CD2BP2 (N Terminus)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 50kDa band observed in Human Lymph Node, Thymus and Tonsil lysates (calculated MW of 37.6kDa according to NP_006101.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Nishizawa et al, Proc Natl Acad Sci U S A. 1998 Dec 8;95(25):14897-902). Recommended concentration: 0.3-1.0 µg/ml.
Protein Sequence	PKRKVTFQGVGDE
Clonality	Polyclonal
MW	37.6
Entrez	10421
Dilution Range	ELISA: 1:32000, WB: 0.3-1.0 µg/ml, IHC-P: 5-10 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Thymus ly...



Immunohistochemical staining of Human Sp...