

Product Datasheet

**phosphor-PPP3CA (phospho-Ser469) antibody
(FITC) [Out of stock] (orb190202)**

Description

Rabbit polyclonal antibody

(Ser469) conjugated to FITC

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

FITC

Tested

ICC

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human PPP3CA around the phosphorylation site of Ser469 (421-521/521aa)

Preservatives

0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

57kDa

Uniprot ID
Q08209
Dilution Range

ICC: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.