

Product Datasheet

FFAR1 antibody (orb19022)

Description	Goat polyclonal antibody to FFAR1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GPR40
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 28kDa band observed in Human Breast lysates (calculated MW of 31.5kDa according to NP_005294.1). Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	CAARTQGGSQK
Clonality	Polyclonal
MW	31.5
Entrez	2864
Dilution Range	ELISA: 1:16000, WB: 0.01-0.03 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Breast ly...