


Product Datasheet

PRDM11 antibody (orb19024)

Description	Goat polyclonal antibody to PRDM11
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	PRDM11
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 50kDa band observed in Human Lung Carcinoma A549 cell line lysates (calculated MW of 58kDa according to NP_064614.2). Recommended concentration: 2-3 µg/ml.
Protein Sequence	LSEGRVRSGLCGG
Clonality	Polyclonal
MW	57.9
Entrez	56981
Dilution Range	ELISA: 1:32000, WB: 2-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb19024 (2ug/ml)
staining of A549
lysat...



Immunohistochemical
staining of Human
Co...