



Product Datasheet

SERPINI1 antibody (orb19025)

Description	Goat polyclonal antibody to SERPINI1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Target	Neuroserpin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 48kDa band observed in Human Brain lysate (calculated MW of 46.4kDa according to NP_005016.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	TMNTSGHDFEEL
Clonality	Polyclonal
MW	46.4
Entrez	5274
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb19025 (1µg/ml)
staining of Human
Brai...



Immunohistochemical
staining of Human
Co...