

Product Datasheet

SERPINI2 antibody (orb19038)

Description

Goat polyclonal antibody to SERPINI2

Species/Host

Goat

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Target

MEPI / SERPINI2

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.IHC: In paraffin embedded Human Pancreas, shows cytoplasmic staining of acinar cells and of rare cells in the islets of Langerhans. Recommended concentration, 3µg/ml.WB: Approx 40kDa band observed in Human Pancreas lysates (calculated MW of 46.1kDa according to NP_006208.1 and 47.4kDa according to NP_001012303.1). Recommended concentration: 0.3-1 µg/ml.

Protein Sequence

KGDWKQKFRKEDTQ

Clonality

Polyclonal

MW

46.1; 47.4

Entrez

5276

Dilution Range

ELISA: 1:64000, WB: 0.3-1 µg/ml, IHC-P: 3µg/ml


Expiration Date

12 months from date of receipt.

250kDa
150kDa
100kDa

50kDa
37kDa
25kDa
20kDa
15kDa
10kDa

Western blot analysis
of Human Pancreas ...



Immunohistochemical
staining of Human
Pa...