


Product Datasheet

HNF4A antibody (orb19044)

Description	Goat polyclonal antibody to HNF4A
Species/Host	Goat
Reactivity	Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	HNF4A
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 52kDa and 48kDa bands observed in lysates of Human hepatocellular carcinoma cell line HepG2 (calculated MW of 52.8kDa according to NP_000448.3 and 46.6kDa according to NP_849181.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	RLSKTLVDMDMADY
Clonality	Polyclonal
MW	51.6; 52.8; 46.6
Entrez	3172
Dilution Range	ELISA: 1:32000, WB: 0.1-0.5 µg/ml, IHC-P: 2 µg/ml
Expiration Date	12 months from date of receipt.



orb19044 (0.1ug/ml)
staining of HepG2
ly...



Immunohistochemical
staining of Human
Sm...



Western blot analysis
of Mouse Small Int...