


Product Datasheet

KLRK1 antibody (orb19046)

Description	Goat polyclonal antibody to KLRK1
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Target	NKG2D / KLRK1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 35kDa band observed in human lymph node and human spleen lysates (calculated MW of 25.3kDa according to NP_031386.1). The observed molecular weight is explained by glycosylation (Han et al, Blood. 2004 Nov 1;104(9):2858-66. Epub 2004 Jul 6; PMID: 15238421). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	KVYSKEDQDLLK
Clonality	Polyclonal
MW	25.3
Entrez	22914
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Spleen ly...