



Product Datasheet

CFH antibody (orb19049)

Description	Goat polyclonal antibody to CFH
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Target	Complement factor H
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 170kDa band observed in human lung lysates (calculated MW of 139kDa according to NP_000177.2). Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	HLVPDRKKDQYK
Clonality	Polyclonal
MW	139
Entrez	3075
Dilution Range	ELISA: 1:32000, WB: 0.03-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb19049 staining (0.03ug/ml) of human l...



orb19049 (5µg/ml) staining of paraffin ...