

Product Datasheet

SP1 antibody (orb19050)

Description	Goat polyclonal antibody to SP1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 100kDa band observed in Human HeLa lysates (calculated MW of 80.7kDa according to NP_612482.2). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	CSRIESPNENSNSQ
Clonality	Polyclonal
MW	80.7; 79.9
Entrez	6667
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb19050
staining
(0.3ug/ml) of
HeLa lys...