

Product Datasheet

SART3 antibody (orb19051)

Description	Goat polyclonal antibody to SART3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SART3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 170kDa band observed in lysates of cell lines Daudi, Jurkat and MOLT4 (calculated MW of 111kDa according to NP_055521.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Yang et al, Cancer Res. 1999 Aug 15;59(16):4056-63; PMID: 10463607). Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	MSNADFAKLFLRK
Clonality	Polyclonal
MW	110
Entrez	9733
Dilution Range	ELISA: 1:16000, WB: 0.01-0.03 µg/ml
Expiration Date	12 months from date of receipt.



orb19051
staining
(0.01ug/ml)
of Jurkat ...