

Product Datasheet

RAB11FIP1 antibody (orb19056)

Description

Goat polyclonal antibody to RAB11FIP1

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Target

RAB11FIP1 / RCP

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:64000. WB: Preliminary experiments gave an approx 90kDa band in lysates of cell line HeLa after 0.1 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 71.0kDa according to NP_079427.3. The 90kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

EATKEAKESKKPE

Clonality

Polyclonal

MW

71.0; 137

Entrez

80223

Dilution Range

ELISA: 1:64000, WB: 0.1-0.3 µg/ml

Expiration Date

12 months from date of receipt.

250kDa
150kDa
100kDa

50kDa
37kDa
25kDa
20kDa
15kDa

Western blot
analysis of
HeLa lysate
usi...