



Product Datasheet

ERLIN2 antibody (orb19057)

Description	Goat polyclonal antibody to ERLIN2
Species/Host	Goat
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	SPFH2 / ERLIN2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 43kDa band observed in Human Breast lysates and in lysates of Human Breast cancer (calculated MW of 37.8kDa according to NP_009106.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Pearce et al, J Biol Chem. 2007 Jul 13;282(28):20104-15.; PMID: 17502376). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	EDEPLETATKEN
Clonality	Polyclonal
MW	37.8
Entrez	11160
Dilution Range	ELISA: 1:2000, WB: 0.3-1 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



orb19057 staining
(1ug/ml) of Human
Brea...



Immunohistochemical
staining of Human
Br...