

Product Datasheet

Recombinant Human Apolipoprotein-Serum Amyloid A (orb1905940)

Description	Recombinant Human Apolipoprotein-Serum Amyloid A
Form/Appearance	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris-HCl, pH 9.0, 150 mM NaCl.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- Refer to lot specific COA for the Use by Date when stored at -20 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Note	For research use only
Protein Sequence	RSFFSFLGEA FDGARDMWRA YSDMREANYI GSDKYFHARG NYDAAKRGPG GVWAAEAISN ARENIQRFFG RGAEDSLADQ AANEWGRSGK DPNHFRPAGL PEKY
Purity	> 98 % by SDS-PAGE and HPLC analyses.
Source	Escherichia coli
MW	Approximately 11.7 kDa, a single non-glycosylated polypeptide chain containing 104 amino acids.
Uniprot ID	P0DJ18
Expiration Date	6 months from date of receipt.

