

Product Datasheet

Recombinant Human Neuroserpin (orb1905969)

Description	Recombinant Human Neuroserpin
Endotoxins	Less than 1 EU/ug of rHuNeuroserpin as determined by LAL method.
Form/Appearance	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Note	For research use only
Protein Sequence	TGATFPEEAI ADLSVNMYNR LRATGEDENI LFSPLSIALA MGMMELGAQG STQKEIRHSM GYDSLKNGEE FSFLKEFSNM VTAKESQYVM KIANSLFVQN GFHVNEEFLQ MMKKYFNAAV NHVDFSQNVV VANYINKWVE NNTNNLVKDL VSPRDFDAAT YLALINAVYF KGNWKSQFRP ENTRTFSFTK DDESEVQIPM MYQQGEFYVG EFS DGSNEAG GIYQVLEIPY EGDEISMMLV LSRQEVPLAT LEPLVKAQLV EEWANSVKKQ KVEVYLPRFT VEQEIDLKDV LKALGITEIF IKDANLTGLS DNKEIFLSKA IHKSFLEVNE EGSEAAVSG MIAISRMVAVL YPQVVDHPF FFLIRNRRTG TILFMGRVMH PETMNTSGHD FEEL
Purity	> 95 % by SDS-PAGE and HPLC analyses.
Source	Escherichia coli
MW	Approximately 44.7 kDa, a single non-glycosylated polypeptide chain containing 394 amino acid residues.
Uniprot ID	Q99574
Expiration Date	6 months from date of receipt.

