


Product Datasheet

**Recombinant Murine Regulation upon Activation
Normal T cell Express Sequence/CCL5
(orb1906023)**

Description

Recombinant Murine Regulation



Endotoxins

Express Sequence/CCL5
Less than 1 EU/ug of rMuRANTES/CCL5 as determined by LAL method.

Form/Appearance

Lyophilized from a 0.2 µm filtered concentrated solution in 30 % Acetonitrile and 0.1 % TFA.

Storage

Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.

Note

For research use only

Protein Sequence

SPYGSDTTPC CFAYLSLALP
RAHVKEYFYT SSKCSNLAVV
FVTRRRQVC ANPEKKWVQE
YINYLEMS

Purity

> 97 % by SDS-PAGE and HPLC analyses.

Source

Escherichia coli

MW

Approximately 7.9 kDa, a single non-glycosylated polypeptide chain containing 68 amino acids.

Uniprot ID

P30882

Expiration Date

6 months from date of receipt.