

Product Datasheet

Recombinant Murine I-TAC/CXCL11 (orb1906033)

Description	Recombinant Murine I-TAC/CXCL11
Endotoxins	Less than 1 EU/ug of rMul-TAC/CXCL11 as determined by LAL method.
Form/Appearance	Lyophilized from a 0.2 µm filtered concentrated solution in 10 mM Sodium Citrate, pH 4.0, with 600 mM NaCl.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Note	For research use only
Protein Sequence	FLMFKQGRCL CIGPGMKAVK MAEIEKASVI YPSNGCDKVE VIVTMKAHKR QRCLDPRSKQ ARLIMQAIEK KNFLRRQNM
Purity	> 98 % by SDS-PAGE and HPLC analyses.
Source	Escherichia coli
MW	Approximately 9.1 kDa, a single non-glycosylated polypeptide chain containing 79 amino acids.
Uniprot ID	Q9JHH5
Expiration Date	6 months from date of receipt.