

Product Datasheet

Recombinant Human Mucosae-associated Epithelial Chemokine/CCL28 (orb1906047)

Description	Recombinant Human Mucosae-associated Epithelial Chemokine/CCL28
Endotoxins	Less than 1 EU/ug of rHuMEC/CCL28 as determined by LAL method.
Form/Appearance	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM PB, pH 7.4, 130 mM NaCl.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Note	For research use only
Protein Sequence	SEAILPIASS CTEVSHHIS RRLLEVNMC RIQRADGDCD LAAVILHVKR RRICVSPHNH TVKQWMKVQA AKKNGKGNVC HRKKHHGKRN SNRAHQGKHE TYGHKTPY
Purity	> 97 % by SDS-PAGE and HPLC analyses.
Source	Escherichia coli
MW	Approximately 12.4 kDa, a single non-glycosylated polypeptide chain containing 108 amino acids.
Uniprot ID	Q9NRJ3
Expiration Date	6 months from date of receipt.

