

Product Datasheet

Recombinant Human Cuteaneous T-cell Attracting Chemokine/CCL27 (orb1906048)

Description
Recombinant Human Ccutaneous

	Chemokine/CCL27
Endotoxins	Less than 0.1 EU/ug of rHuCTACK/CCL27 as determined by LAL method.
Form/Appearance	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Note	For research use only
Protein Sequence	FLLPPSTACC TQLYRKPLSD KLLRKVIQVE LQEADGDCHL QAFVLHLAQR SICIHPQNPS LSQWFEHQER KLHGTLPLKLN FGMLRKMG
Purity	> 96 % by SDS-PAGE and HPLC analyses.
Source	Escherichia coli
MW	Approximately 10.1 kDa, a single, non-glycosylated polypeptide chain containing 88 amino acids.
Uniprot ID	Q9Y4X3
Expiration Date	6 months from date of receipt.