

Product Datasheet

Recombinant Human LD78 beta/CCL3L1 (orb1906071)

Description	Recombinant Human LD78 beta/CCL3L1
Endotoxins	Less than 1 EU/ug of rHuLD78-β/CCL3L1 as determined by LAL method.
Form/Appearance	Lyophilized from a 0.2 μm filtered concentrated solution in 20 mM PB, pH 7.4, 100 mM NaCl.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Note	For research use only
Protein Sequence	APLAADTPTA CCFSYTSRQI PQNFIADYFE TSSQCSKPSV IFLTKRGRQV CADPSEEWVQ KYVSDLELSA
Purity	> 97 % by SDS-PAGE and HPLC analyses.
Source	Escherichia coli
MW	Approximately 7.8 kDa protein containing 70 amino acid residues, including the four highly conserved cysteine residues present in CC chemokines.
Uniprot ID	P16619
Expiration Date	6 months from date of receipt.

