

Product Datasheet

Recombinant Porcine Interleukin-8/CXCL8 (orb1906096)

Description	Recombinant Porcine Interleukin-8/CXCL8
Endotoxins	Less than 1 EU/ug of rPoIL-8/CXCL8 as determined by LAL method.
Form/Appearance	Lyophilized from a 0.2 µm filtered concentrated solution in 2 x PBS, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Note	For research use only
Protein Sequence	ARVSAELRCQ CINTHSTPFH PKFIKELRVI ESGPHCENSE IIVKLVNGKE VCLDPKEKWV QKVVQIFLKR TEKQQQQQ
Purity	> 95 % by SDS-PAGE and HPLC analyses.
Source	Escherichia coli
MW	Approximately 9.1 kDa, a single non-glycosylated polypeptide chain containing 78 amino acids.
Uniprot ID	P26894
Expiration Date	6 months from date of receipt.

