


Product Datasheet

SDHB antibody (orb19061)

Description	Goat polyclonal antibody to SDHB
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	SDHB
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 28kDa band observed in Human Liver lysates (calculated MW of 31.6kDa according to NP_002991.1). Recommended concentration: 0.03-0.1 µg/ml.
Protein Sequence	ATYKEKKASV
Clonality	Polyclonal
MW	31.6
Entrez	6390
Dilution Range	ELISA: 1:8000, WB: 0.03-0.1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Liver lys...



Immunohistochemical staining of Human He...