

Product Datasheet

Recombinant Human Heparin-binding EGF-like Growth Factor (orb1906266)

Description	Recombinant Human Heparin-binding EGF-like Growth Factor
Endotoxins	Less than 0.1 EU/ug of rHuHB-EGF as determined by LAL method.
Form/Appearance	Lyophilized from a 0.2 µm filtered concentrated solution in 10 mM PB, 500 mM NaCl, 5% trehalose, pH 7.4.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Note	For research use only
Purity	> 97% by SDS-PAGE and HPLC analyses.
Source	Escherichia coli
MW	Approximately 9.7 kDa, a single non-glycosylated polypeptide chain containing 86 amino acids.
Uniprot ID	Q99075
Expiration Date	6 months from date of receipt.

