

Product Datasheet

Cytokeratin 20 Antibody (orb1906367)

Description

Cytokeratin 20 Antibody

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The original antibody was generated by immunizing female ICR mice with human CK20 peptide (amino acids 411-424) conjugated to KLH in equal volume of Freund's incomplete adjuvant.
Target	Cytokeratin 20
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody can bind human cytokeratin 20 in an ELISA. It is also capable of detecting human cytokeratin 20 in a western blot assay. This antibody is similar to commercially available anti-Cytokeratin 20 antibody KS20.8 in staining various human tissues. This antibody was used in the immunohistochemical staining of paraffin sections of stomach, colon and colonic cancer tissues (CN109734805).
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P35900
Expiration Date	12 months from date of receipt.