

Product Datasheet

Recombinant SARS-CoV-2 3C-like Proteinase (C-6His) (orb1906479)

Description	Recombinant SARS-CoV-2 3C-like Proteinase is produced by our E.coli expression system and the target...
Form/Appearance	Supplied as a 0.2 µm filtered solution of 50mM Tris-HCl, 1mM EDTA, pH 7.3.
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-6His
Note	For research use only
Protein Sequence	Ser3264-Gln3569
Source	E.coli
Uniprot ID	P0DTC1
Expiration Date	6 months from date of receipt.