

## Product Datasheet

### Recombinant Human B7-1/CD80 (C-6His) (orb1906493)

|                         |                                                                                                                                                                                                                                          |
|-------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>      | Recombinant Human Activation B7-1 antigen is produced by our Mammalian expression system and the tar...                                                                                                                                  |
| <b>Form/Appearance</b>  | Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.                                                                                                                                                                              |
| <b>Storage</b>          | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks.<br>Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| <b>Tag</b>              | C-His                                                                                                                                                                                                                                    |
| <b>Note</b>             | For research use only                                                                                                                                                                                                                    |
| <b>Protein Sequence</b> | Val35-Asn242                                                                                                                                                                                                                             |
| <b>Source</b>           | Human Cells                                                                                                                                                                                                                              |
| <b>Uniprot ID</b>       | <b>P33681</b>                                                                                                                                                                                                                            |
| <b>Expiration Date</b>  | 6 months from date of receipt.                                                                                                                                                                                                           |